

```

# freq_modele.py

01| import csv
02| import os
03| import inspect
04|
05| from detecta import detect_peaks
06|
07| ##Débit moyen
08| file=open('C:\\Users\\Louis\\Documents\\TIPE\\modélisation2\\5Lmin\\signal.txt',
09| 'r')
10| reader = csv.reader(file,delimiter='\\t',skipinitialspace=True)
11|
11| X = []
12| Y = []
13| A = []
14|
15| for row in reader:
16|     X.append(float(row[0]))
17|     Y.append(float(row[1]))
18|
19| moy = 0
20|
21| for u in range(len(Y)-1):
22|     A.append((Y[u]+Y[u+1])/2*(X[u+1]-X[u]))
23|
24| for p in range(len(A)):
25|     moy = moy + A[p]
26| moy = moy/1000
27|
28| ##max min
29| file2=open('C:\\Users\\Louis\\Documents\\TIPE\\modélisation2\\5Lmin\\
30| \\pression_model.txt', 'r')
31| reader = csv.reader(file2,delimiter=',',skipinitialspace=True)
32|
32| XP = []
33| YP = []
34|
35| for row in reader:
36|     XP.append(float(row[0]))
37|     YP.append(float(row[1]))
38|
39| maxP = max(YP)
40| minP = min(YP[XP.index(5):]) #va chercher la valeur min dans YP à partir de la
41| valeur XP[i]=5 jusqu'à la fin de la liste
42|
42| DP = maxP-minP
43|
44| ##détection fréquence cardiaque
45| i = XP.index(int(3))
46| YPFC = YP[i:]
47| XPFC = XP[i:]
48|
49| pic = detect_peaks(YPFC)
50|
51| picok = []
52| m = len(pic)
53|
54| for l in range(m):
55|     if pic[l]>maxP - 5:
56|         picok.append(pic[l])
57|
58|
59| DT = XPFC[picok[-1]]-XPFC[picok[0]]
60|
61| FC = (len(picok)-1)*60/DT
62|

```

```
63| print("Débit cardiaque:", round(moy,3), "L/min \n Delta pression", round(DP, 3),  
"mmHg \n Pression max", round(maxP, 3) , "mmHg \n Pression min", round(minP, 3), "mmHg  
\n Fréquence cardiaque", round(FC, 3), "bat/min")
```